



Product Description: HoMedics Thera-P Model HP-200 Standard Size Moist/Dry Heating Pad with Auto Shut-off and HoMedics Thera-P Model HP-300 King Size Moist/Dry Heating Pad

Each HoMedics heating pad is marked with a unique 4-digit date code located both on the back of the hand control as well as on the bottom panel of the color box. ONLY 4-DIGIT DATE CODES ENDING IN "01" ARE SUBJECT TO THIS RECALL.

In order to accurately identify the date code, consumers should refer to the diagram included in this release.

Manufacturer: HoMedics, Inc., Commerce Township, Michigan

Date Recall Posted: 2/21/07

UPCs: Meijer has not carried this product since 2001. UPC Numbers currently not available in the Meijer system.

Summary Of Recall: It has been determined that some of the heating pads contained an inadequate connector crimp, which lead to a high resistance connection that generated excessive heat, thereby posing a risk of burn injuries, fire or damage to the heating pad itself or to materials (like bedding and furniture) that could come into contact with the pad.

What You Should Do: Consumers whom have any of these particular heating pads should discontinue the use of the product immediately and return it to the retailer of purchase for a full refund. HoMedics is currently contacting the retailers of these heating pads to communicate return instructions. Further information can be obtained by going to www.homedics.com.

Consumers with questions may contact the company at 1-800-466-3342.